

HIGH RESOLUTION PULSE OXIMETRY(HRPO) CASE REPORT

Hypoxemia/
Hypoventilation
Syndrome Example

PHYSICIAN:

DOB:

HEIGHT (in): -

BMI: -

WEIGHT (lb): -

NECK SIZE:


UPLOADED AUTOMATED FROM PROFOX FILE

This indicator examines the overall trend of oxygen levels independent of events, and highlights gradual drifts to lower baseline oxygen saturation levels.

START:

END:

DURATION:



PATTERN BASED REPORT

OXYGEN SATURATION BASELINE ANALYSIS


Oxygen Saturation Baseline Drift(OSBG) (normal <= 3)	15
Initial Saturation Baseline	93
Lowest Saturation Baseline	81
Highest Saturation Baseline	97

Baseline is determined by the Mean SpO2 during 3 Minute window without Artifact and without Events.

Period SpO2 was < 88%

The total time oxygen saturation was <= 88% was 00:13:36

DESAT SEVERITY
(most severe 10%)



SPO2 CYCLING

% Time in Cycling (Duration)	8%	(00:37:36)
Cycling Frequency	22	
96% - Lowest Sat	10	
Cycling Severity Index	9	

Pattern based report gives detailed information on desaturations & cyclical patterns.

TRADITIONAL REPORT

ODI4:	3
Total OD4 Events:	27
Time in OD4 Events:	00:40:30
Avg OD4 Event Duration:	00:00:51
<=88% OD4 Events:	3
<=88% Longest Duration:	00:10:09
Minimum SpO2:	80
Avg Low 10% SpO2:	84
Avg Low SpO2:	91
Avg Low SpO2 <=88%:	84

Definition of OD4 Event: a fall in oxygen saturation of at least 4% and persisting greater than 3 seconds.

%SpO2	DURATION	%TOTAL
94-100	05:19:12	64%
88-94	02:47:29	33%
80-88	00:13:34	3%
70-80	00:00:02	0%
<= 70	00:00:00	0%
Total	08:20:17	100%
Motion Artifact	00:00:19	0.06%
Error Signal	00:00:01	0%

Traditional oximetry report parameters.

HIGH RESOLUTION OXIMETRY CASE REPORT





PATIENT ID:

START:


END:

Subsequent pages display event examples and raw collection data.


Example of Detected Cycling




Percent of study spent in the specified SpO2 range.




Event count for the specified SpO2 range.




MOTION ARTIFACT
ERROR SIGNAL



04:10 02:30 12:50 04:10



MOTION ARTIFACT
ERROR SIGNAL



04:10 05:50 07:30 09:10


Detailed raw data report.
(2 hr view)

XIMETRY CASE REPORT


PATIENT ID:

PAGE START TIME:


PAGE DURATION: 2 hours




MOTION ARTIFACT
ERROR SIGNAL



MOTION ARTIFACT
ERROR SIGNAL



MOTION ARTIFACT
ERROR SIGNAL



MOTION ARTIFACT
ERROR SIGNAL